# Tech-Notes

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
Figure 1. Annealed microstructure of type O6 graphitic tool steel prepared (top) with the three-step method using the UltraPol® silk cloth for step 2; and (bottom) using the ApexHercules® H rigid grinding disc for step 2 (magnification bars are 10-μm long; 4% picral etch).

Table 1. Four-step method for preparing tool steels

Abrasive & Surface	Lubricant	RPM	Head Platen/ Direction	Load per Specimen	Time (minutes)
120- to 240-grit* (P120 to P280) SiC CarbiMet® 2 waterproof abrasive paper, or 125 to 45μmApex® DGDs or DGD Red or Purple	water	240-300	Contra	6 lbs (27N)	Until Plane
9-μm diamond on an UltraPol® silk cloth (or UltraPad® polyester cloth or ApexHercules® H rigid grinding disc)	MetaDi® Fluid	120-150	Contra	6 lbs (27N)	5
3-μm diamond on a TriDent® cloth (or TexMet® chemotextile pad)	MetaDi® Fluid	120-150	Contra	6 lbs (27N)	3
MasterPrep® Alumina Suspension on a MicroCloth® pad	No other lubricant is needed	120-150	Contra	6 lbs (27N)	1-3

\* Use 120-grit for specimens > 60 HRC; use 180-grit for specimens at 35-60 HRC; use 240-grit for specimens < 35 HRC.

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
Figure 2. Spheroidize annealed W1 water-hardened carbon tool steel etched with (top) 4% picral to reveal the general structure, with (middle) Klemm's I to color the ferritic matrix, and with (bottom) alkaline sodium picrate at 90 °C for 60 seconds to color the cementite.

Table 2. Three-step method for preparing tool steels

Abrasive & Surface	Lubricant	Platen RPM	Head Platen/ Direction	Load per Specimen	Time (minutes)
120- to 240-grit (P120 to P280) SiC CarbiMert® 2 waterproof abrasive paper* or 125 to 45μm Apex® DGDs or DGD Red or Purple	water	240-300	Contra	6 lbs (27N)	Until Plane
3-μm diamond on an UltraPol® silk cloth cloth (or UltraPad® polyester cloth or ApexHercules® H rigid grinding disc)	MetaDi® Fluid	120-150	Contra	6 lbs (27N)	5
MasterPrep® Alumina Suspension on a MicroCloth® pad	No other lubricant is needed	120-150	Contra	6 lbs (27N)	5

\* Use 120-grit for specimens > 60 HRC; use 180-grit for specimens at 35-60 HRC; use 240-grit for specimens < 35 HRC.

Table 3. Etchants for Tool Steels

Composition	Comments
1-10 mL HNO <sub>3</sub> 99-90 mL Ethanol	Nital. Most commonly used reagent. Reveals ferrite grain boundaries and ferrite-carbide interfaces. Excellent for martensite. Do not store solutions with >3% HNO <sub>3</sub> in a tightly stopped bottle.
4 g Picric acid 100 mL Ethanol	Picral. Recommended for annealed structures or those containing pearlite or bainite. Does not reveal ferrite grain boundaries. Addition of a few drops of zephiran chloride increases etch rate. Add 1-5 mL HCl to improve etch response for annealed higher alloy tool steels.
1 g Picric acid 5 mL HCl 95 mL Ethanol	Vilella's reagent. Reveals structure of higher alloyed tool steels.
50 mL sat. Aqueous Sodium thiosulfate 1 g Potassium metabisulfite	Klemm's I tint etch. Immerse specimen until the surface is colored violet. Colors ferrite blue and red while martensite is brown. Carbides are unaffected. Works well only on low alloy and carbon tool steels.
0.6 mL HCl 0.5-1.0 g Potassium metabisulfite 100 mL water	Beraha's reagent. Immerse specimen until the surface is colored. Colors ferrite and martensite, carbides are not affected. Good for most tool steels.
3 g Potassium metabisulfite 2 g Sulfamic acid 0.5-1.0 g Ammonium bifluoride 100 mL water	Beraha's sulfamic acid reagent No. 4. For carbon and low-alloy tool steels, leave out the NH <sub>4</sub> F·HF. Immerse until the surface is colored. Ferrite and martensite are colored; carbides are not affected. Good for high chromium tool steels.
2 g Picric acid 25 g NaOH 100 mL water	Alkaline sodium picrate. Colors cementite and M <sub>6</sub> C carbides. Immerse specimen in solution at 80-100 °C for 1 minute or more.
10 g K <sub>3</sub> Fe(CN) <sub>6</sub> 10 g NaOH or KOH 100 mL water	Murakami's reagent. Use at 20 °C to outline and darken M <sub>7</sub> C <sub>3</sub> and M <sub>6</sub> C, and to outline M <sub>2</sub> C. M <sub>23</sub> C <sub>6</sub> is faintly colored.
4 g KMnO <sub>4</sub> 4 g NaOH 100 mL water	Groesbeck's reagent. Use at 20 °C to outline M <sub>2</sub> C and to outline and darken M <sub>6</sub> C. M <sub>7</sub> C <sub>3</sub> is faintly colored.
1 g CrO <sub>3</sub> 100 mL water	Blickwede and Cohen's etch. Use at 2-3 V dc, 20 °C, for 30 seconds with a stainless steel cathode. Outlines M <sub>23</sub> C <sub>6</sub> , outlines and colors M <sub>7</sub> C <sub>3</sub> , colors MC and attacks M <sub>2</sub> C.

Table 4. Results of the Etching Experiments

Etchant	M <sub>3</sub> C	M <sub>23</sub> C <sub>6</sub>	M <sub>7</sub> C <sub>3</sub>	M <sub>6</sub> C	MC	M <sub>2</sub> C
Alk.Na Pic.	Colors	NA	NA	Colors	NA	NA
Murakami	NA	Faint	Outlined/ Colored	Outlined/ Colored	NA	Outlined
Groesbeck	NA	NA	Faint	Outlined/ Colored	NA	Outlined
1% CrO <sub>3</sub>	NA	Outlined	Outlined/ Colored	NA	Colors	Attacks
NA – no affect						

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Thomas Li  
[Thomas.li@buehler.com](mailto:Thomas.li@buehler.com)  
Room C0205, Building C, College of materials science and engineering, Tianjing University, China  
Web Site: <http://www.buehler-asia.com>

**BUEHLER**  
The Science Behind Materials  
Preparation & Analysis

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BUEHLER Ltd.